



Data sheet Monkey IL-1 β ELISPOT antibody pair

Cat. No.: ACT607-20

Coating antibodies (1 vial)

Product: Monoclonal antibody to monkey interleukin 1 beta (IL-1 β)
Isotype: Mouse IgG₁
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Contents: Each vial contains sufficient material for coating of twenty 96-well ELISPOT plates
Buffer: Prior to lyophilization: 1.0 ml PBS + 125 mM trehalose
Application: Coating antibody in an ELISPOT system
Reconstitution: Dissolve the contents of the vial by injection of 1.0 ml distilled water into the vial and dilute 100 times in PBS. The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 50 μ l/well).

Detection antibodies (1 vial)

Product: Biotinylated monoclonal antibody to monkey interleukin 1 beta (IL-1 β)
Isotype: Mouse IgG₁
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)
Contents: Each vial contains sufficient material for twenty 96-well ELISPOT plates
Buffer: Prior to lyophilization: 2.0 ml PBS + 1% BSA + 125 mM trehalose
Application: Detection antibody in an ELISPOT system
Reconstitution: Dissolve the contents of the vial by injection of 2.0 ml distilled water into the vial and dilute 100 times in Dilution buffer (see TDS). The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 100 μ l/well).

General

Specificity: Validated for detecting rhesus macaque, pig-tailed macaque and lion-tailed macaque IL-1 β
Sterility: Membrane filtered (0.2 μ m)
Stability: The lyophilized products are stable for more than one year at 4°C.
After reconstitution, the antibodies are stable for several months at 4°C (if kept sterile) or for minimal one year at -20°C.