



Manufactured By: HYPHEN BioMed

Anti-Rat-PAI-1 Coupled with Peroxidase
Anti-Rat-PAI-1 peroxydase
Vial of about 16 µg/flacon d'environ 16 µg

Ref. AAP009A

STORE AT
2-8°C

For research use only

Uniquement à usage de recherche

CONSERVER A
2-8°C

Lot :

Exp :

Presentation:

Vial containing about 16 µg of immunoconjugate, anti-rat-PAI-1 Peroxidase, lyophilised.

Reconstitution:

Restore the vial with 7.5 ml of ELISA dilution buffer.

Immunoconjugate:

Monoclonal antibody specific for rat-PAI-1 conjugated to peroxidase.

Use:

For in vitro and research use only.

Second antibody for the revelation of insolubilised or immunocapture rat-PAI-1

Stability:

Restored : - 7 days at 2-8°C

Lyophilized : - Until the expiration date printed on the vial.

Présentation :

Flacon lyophilisé contenant environ 16 µg d'immunoconjugué, anti-PAI-1 de rat peroxydase.

Reconstitution:

Chaque flacon doit être reconstitué par 7.5 ml de tampon de dilution ELISA

Immunoconjugué :

Anticorps monoclonal anti-PAI-1 de rat couplé à la peroxydase.

Utilisation :

Uniquement à usage de recherche, in vitro.

Second anticorps permettant de révéler le PAI-1 de rat insolubilisé, ou immobilisé.

Stabilité du flacon:

Reconstitué : - 7 jours à 2-8°C.

Lyophilisé : - Jusqu'à la date d'expiration imprimée sur le flacon.

Rev 080402
D.750.08/AP/009A



8560 Gove Court • Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com