



Manufactured By: HYPHEN BioMed

**Anti-TAFI-Peroxidase conjugate**  
**Conjugué anti-TAFI couplé à la peroxydase**  
Vial of 0.5 mg/flacon de 0.5 mg

Ref. AAP011Z

STORE AT  
2-8°C

For research use only



Uniquement à usage de recherche

CONSERVER A  
2-8°C

Lot :

Expiration :

**Origin:**

Highly purified goat antibodies specific for human TAFI coupled to Horse Radish Peroxidase (HRP).

**Presentation:**

Vial containing 0.5 mg of immunoconjugate (1 ml at 500 µg/ml), in the liquid form.

**Excipients:**

Bovine Serum Albumin (BSA), Sodium chloride, phosphate buffer, stabilizers. Kathon CG is present at a concentration of 0.05%.

**Use:**

For *in vitro* use only.

Can be used as a tag antibody in ELISA, immuno-precipitation, immuno-cytochemistry or immunoblotting, for studies on TAFI.

The usual working dilutions are in the range 1:100 to 1:1000, in a dilution buffer containing a carrier protein (Bovine Serum Albumin or casein for example).

**Stability:**

In its original vial, until the expiration date printed on the vial.

Store at 2-8°C.

Avoid any contamination when opening the vial.

**Origine :**

Anticorps de chèvre anti-(h) TAFI hautement purifié et couplé à la peroxydase (HRP).

**Présentation :**

Flacon contenant 0.5 mg d'immunoconjugué (1 ml at 500 µg/ml), sous forme liquide.

**Excipients :**

Albumine bovine (BSA), Chlorure de sodium, tampon phosphate, stabilisants. Le réactif contient 0,05% de Kathon CG.

**Utilisation :**

Utilisation *in vitro* exclusivement.

Peut être utilisé comme anticorps marqué en ELISA, immuno-précipitation, immuno-cytochimie ou immunoblotting, pour les études sur TAFI.

La dilution habituellement utilisée est comprise entre 1/100 et 1/1000, dans un diluant contenant un ballast protéique (albumine bovine ou caséine par exemple).

**Stabilité du flacon:**

Conservé dans son flacon d'origine, le produit est stable jusqu'à la date d'expiration imprimée sur le flacon.

Conserver à 2-8°C

Eviter toute contamination à l'ouverture du flacon.

AP011Z/3  
F.750.08/AAP/011Z



6560 Gove Court • Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com