



Manufactured By: HYPHEN BioMed

**ELISA Conjugate Diluent buffer, with albumin  
Tampon de dilution pour Immunoconjugué  
(albuminé)  
(Vial of 25 ml/flacon de 25 ml)**

Ref. AAR017A

STORE AT  
2-8°C

For research use only

Lot:



BioHazard

Uniquement à usage de recherche

Expiration:

CONSERVER A  
2-8°C

**Composition:**

Phosphate buffer saline (9 g/l NaCl) containing bovine serum albumin (1%) and Tween 20 (0.1%).

Contains Kathon CG as preservative.

Ready to use.

**pH:**

7.50.

**Use:**

Conjugate Diluent buffer for enzyme immunoassays.

Ready to use

For in vitro use only.

**Stability:**

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

**Composition :**

Tampon phosphate isotonique (9 g/l NaCl) contenant de l'albumine bovine (1%) et du Tween 20 (0.1 %).

Contient du Kathon CG comme conservateur.

Prêt à l'emploi.

**pH :**

7.50

**Utilisation:**

Tampon de dilution des immunoconjugués pour dosages immuno-enzymatiques.

Prêt à l'emploi.

Utilisation *in vitro* exclusivement.

**Stabilité du flacon:**

Dans son flacon d'origine, fermé et conservé à 2-8°C, jusqu'à la date de péremption imprimée sur l'étiquette.

Après ouverture et protégée de toute contamination, cette solution est stable 4 semaines à 2-8°C.

Rev 080328  
D.750.08/AR017A



8580 Gove Court - Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com