



***Cross-Reactivities with Homologous Animal Proteins for the Various ZYMUTEST Kits***

<i>Species/Factor</i>	<i>Human</i>	<i>Monkey Cynomolgus</i>	<i>Monkey Rhésus</i>	<i>Pig</i>	<i>Goat</i>	<i>Sheep</i>	<i>Horse</i>	<i>Dog</i>	<i>Rabbit</i>	<i>Guinea Pig</i>	<i>Rat</i>	<i>Mouse</i>
Annexin V	++++	ND	ND	++++	++++	++++	+++	++	-	+++	+	++++
PF4	++++	++++	++++	+/-	++++	+	-	-	-	+++	++	+++
CRP	++++	++++	++++	-	-	-	-	-	++++	-	-	-
tPA Antigen	++++	++++	++	++++	++++	++++	+++	+	++++	+/-	+/-	+/-
PAI-1 Antigen	++++	++	++++	-	+/-	+/-	-	-	-	-	-	-
Protéine S, Free	++++	++++	++++	-	-	-	-	-	-	-	-	-
Protéine S, Total	++++	++	++	-	-	-	-	-	-	-	-	-
β <sub>2</sub> GPI	++++	++++	++++	+++	+/-	+	++	++	-	+++	+++	++
Fibrinogen	++++	++++	++++	++	+	++	++	+	-	-	+	+
Protein C	++++	++++	++++	+/-	++++	++	-	-	-	+/-	-	-
Fibronectin	++++	+++	++++	++++	++++	++++	++++	++++	-	++++	++++	++++
vWF	++++	++++	++++	++++	++++	++++	++++	++++	-	++	+++	++
Protein Z	++++	++++	++++	+/-	+++	++	+	+	+/-	++	+	++
Factor IX	++++	++++	++++	++++	++++	+++	+++	+++	+/-	++++	++++	++++
Factor X	++++	++++	++++	++	+++	++	+	+	-	+++	++	++
Factor XIII-A	++++	++	++	+++	+++	+	+	+/-	++	-	-	-

**Estimation of respective cross-reactivities:**

++++	> 60 % cross-reactivity
+++	40 to 60 % cross-reactivity
++	20 to 40 % cross-reactivity
+	10 to 20 % cross-reactivity
+/-	2 to 10 % cross-reactivity
-	< 2 % cross-reactivity

Form AH90  
08-2008

**Aniara**

6560 Gove Court  
Mason, OH 45040

**Phone:** 513-770-1991 **Toll Free:** 866-783-3797  
**Web:** www.aniara.com

**FAX:** 513-573-9241  
**Email:** info@aniara.com