

Selection Table for SmartBuffers



8500 Gove Court • Mason, OH 45040
 Phone: 513.770.1901 Toll Free: 866.783.3797 Fax: 513.573.9241 Email: info@anlara.com
www.anlara.com

	PBS	PBS-T	TBS	TBS-T	BBS	NaPi	TAE	TBE	TG	TG+SDS	TE	Carb-Bicarb	SSC	Tris-HCl	NaAc	BSC	EDTA	SSPE	
Applications																			
Nucleic acid electrophoresis							•	•										•	
Protein electrophoresis									•	•									
Western blot		•	•	•	•				•										
Southern blot													•						
Northern blot							•						•						
Coating procedures		•			•							•							
Immunoassays/ ELISA	•	•	•	•	•							•							
Immunohistochemistry	•		•								•								
Microarrays																			•
Tissue and cell culture	•		•											•			•	•	
Microbiology	•																		
Molecular biology			•			•					•		•	•	•				•
Biochemistry/Clinical chemistry	•				•	•								•		•	•	•	
Characteristics																			
Isotonic	•	•	•	•	•								(1x)						
Autoclavable	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Pouches	•		•			•	•	•	•	•	•		•	•	•	•	•	•	
Tablets	•	•	•	•	•							•	•						•
Solution	•																		
Pack size range	16	5	4	2	1	7	2	3	2	2	1	5	3	3	2	2	2	2	1