

ANIARA

Manufactured By: HYPHEN BioMed

Enzyme immunoassays for the exploration of fibrinolysis in human plasma

ZYMUTEST TAFI Ref. ARK008A


Kit presentation

Two site ELISA

Two level quality Controls

Ready to use, complete and disposable kit

Expressed in % of a normal pooled plasma



Assay characteristics

- Total assay time: 3h05 min
- Detection threshold: $\leq 5\%$
- Dynamic range: 1 to ≈ 100 ng/mL
- Intra-assay CV: 3 to 8%
- Inter-assay CV: 5 to 10%

Other products

ZYMUTEST tPA Ag – Ref. ARK011A - CE

ZYMUTEST PAI-1 Ag – Ref. ARK012A - CE

ZYMUTEST PAI-1 Act – Ref. ARK019A

ZYMUTEST uPA – Ref. ARK013A

ZYMUTEST tPA-PAI-1 complexes – Ref. ARK017A

ZYMUTEST uPA-PAI-1 complexes – Ref. ARK018A

Form AH57
5-2007