



Enzyme immunoassays for the exploration of fibrinolysis in human plasma

ZYMUTEST TAFI
Ref. ARK008A


Kit presentation

Two site ELISA

Two level quality Controls

Ready to use, complete and disposable kit

Expressed in % of a normal pooled plasma



Assay characteristics

- Total assay time: 3h05 min
- Detection threshold: $\leq 5\%$
- Dynamic range: 1 to ≈ 100 ng/mL
- Intra-assay CV: 3 to 8%
- Inter-assay CV: 5 to 10%

Other products

ZYMUTEST tPA Ag	Ref. ARK011A - CE
ZYMUTEST PAI-1 Ag	Ref. ARK012A - CE
ZYMUTEST PAI-1 Act	Ref. ARK019A
ZYMUTEST uPA	Ref. ARK013A
ZYMUTEST tPA-PAI-1 complexes	Ref. ARK017A
ZYMUTEST uPA-PAI-1 complexes	Ref. ARK018A

Form AH57
5-2007