

PBS-Tween® tablets

1 tablet in 1000 ml of deionized water yields:

0.010 M Phosphate buffer pH 7.4

0.0027 M KCl

0.140 M NaCl

0.05 % Tween® 20

ANIARA

Manufactured By: Medicago AB

Carbonate-bicarbonate buffer
pH 9.6

One tablet dissolved in 100 ml of water yields
0.05 M Sodium carbonate-bicarbonate buffer pH 9.6

For laboratory use only.

Medicago AB
S-755 98 Uppsala, Sweden
tel +46 18 561180
www.medicago.se



Medicago AB
S-755 98 Uppsala, Sweden
tel +46 18 561180
www.medicago.se

Don't **waste time** preparing your own buffer solutions!

With Medicago's buffer tablets you can focus **on research!**

Medicago's buffer
tablets save time
and money!

304-01

2007-01

+18 to +30

LOT 100501

MADE IN SWEDEN

Medicago AB
S-755 98 Uppsala, Sweden
tel +46 18 561180
www.medicago.se



ANIARA

8580 Gove Court • Mason, OH 45040

Phone: 513.770.1991 Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com

Form AM01
08-2009