

XENOMETRIX


Ames II



Ames MPF™ 98/100

Ames MPF™ 1535

Ames MPF™ 1537



- The Ames MPF and Ames II Fluctuation Assay come as complete kit with strains and liquid culture media ready to use
- The Ames MPF and Ames II Fluctuation Assays can be used to detect genotoxic activity in chemicals, cosmetics, pharmaceuticals, foodstuffs, water, air, soils or sediments.
- Quality controlled *S. typhimurium* strains (genotyped and phenotyped)
- At least 3 fold less test compound as compared to the traditional Ames test.
- Less hands-on-time and 6 times less plastic ware than with the traditional Ames test
- Less consumption of rat liver S9 as compared to the traditional Ames test.
- Microplate format suitable for automation

Form AA19  
03-2009



7768 Service Center Drive • West Chester OH 45069

Phone: 513.770.1991

Toll Free: 866.783.3797


Fax: 513.573.9241

Email: [info@aniara.com](mailto:info@aniara.com)

[www.aniara.com](http://www.aniara.com)

# Ames MPF 98, Ames MPF 100, Ames MPF 1535, Ames MPF 1537, Ames II

## Test procedure



## Please send me more information:

- Instruction Manual:
  Ames MPF 98
  Ames MPF 100
  Ames MPF 1535
  Ames MPF 1537
  Ames II
  Ames MPF 98/100
  Ames MPF 98/100/1535/1537

- Price list  
 References  
 Please contact me

Name: \_\_\_\_\_ First Name: \_\_\_\_\_

Affiliation: \_\_\_\_\_

Street: \_\_\_\_\_

Zip/City/Country: \_\_\_\_\_

Email/Tel.: \_\_\_\_\_