

APROTININ

REF 5D-62106-1

*(Trypsin inhibitor, pancreatic Alkaline)
from bovine lung*

For Laboratory Use Only.
Not for use in human or direct animal applications.

Store at 2-8°C
Lyophilized Powder, 1 g

DESCRIPTION

APROTININ, 1 g, lyophilized powder. Trypsin inhibitor, pancreatic alkaline. Aprotinin is a multiple inhibitor of many serine proteases, including Trypsin, Plasmin, Urokinase, Kallikrein, Factor XIIa, etc. It contains 58 amino acids. It forms stable complexes with these enzymes.

FORMULA

C₂₈₄H₄₃₂N₈₄O₇₉S₇

MOLECULAR WEIGHT

6511

COMPOSITION

Aprotinin powder, lyophilized.

SOLUBILITY

≥ 25 mg/mL in distilled water.

APROTININ ACTIVITY

Aprotinin was tested with Pharmacopoeias' compliant methods for its inhibitory potency of Trypsin, expressed as "Trypsin Inhibitory Units" (TIU). 1 TIU inhibits 50% of the activity of 2 Units of Trypsin.

PRINCIPLE

Aprotinin can be used in any application where inhibition of its targeted proteases is needed.

Current uses include protein isolation, removal of reactive proteases through solid phase, inhibition of undesired proteolytic activity in laboratory assays, or other techniques; for blocking fibrinolytic activity in vitro; etc.

PREPARATION

The vial contains 1 g and can be restored with 10 mL distilled water for obtaining a solution at about 100 mg/mL.

STANDARDIZATION

Aprotinin has a specific activity of about 6,000 KIU/mg (this corresponds to about 6.0 TIU/mg or about 3.35 PEU/mg).

APPLICATIONS

Aprotinin can be used for e.g., the determination of:

1. Targeted proteases inhibition
2. Laboratory assays and techniques
3. Protein purification

STORAGE AND STABILITY

Lyophilized Aprotinin: Stable at 2-8°C up to the expiration date printed on the label; after reconstitution with distilled water: up to 6 months stored at ≤ -20°C.

Aprotinin Solution: contamination by microorganisms must be avoided.

NOTE: Aprotinin contains no preservatives; after reconstitution it must be used rapidly or stored frozen. Alternatively, a preservative (for example Sodium Azide at 0.9 mg/mL) can be added.

WARNINGS AND PRECAUTIONS

For Laboratory Use Only. Not for use in human or direct animal applications.

REFERENCES

1. Lottenberg R, Sjak SN, Fazleabs AT, Roberts RM. Aprotinin inhibits urokinase but not tissue-type plasminogen activator. Thromb Res 1988; 49: 549-56.
2. Fritz H, Wunderer G. Biochemistry and applications of aprotinin, the kallikrein-inhibitor from bovine organs. Arzneimittel Forsch./Drug. Res. 1983; 33: 479-94.



5-Diagnostics AG

Heuberg 7, 4051 Basel,
Switzerland

Tel: +41 61 588 07 84

www.5-diagnostics.com

info@5-diagnostics.com