



IMMUNO REAGENTS

Stop solution for TMB Substrate



Features

- Ready-to-use
- ELISA-reagent
- Sufficient for 10 96-well plates
- Contains 0.5 M H₂SO₄

Product description

Medicago's Stop solution is especially developed for the EC-Blue Enhanced Substrate (TMB) (Article no. 10-9405) and for use in ELISA procedures. The Stop solution is supplied as a ready-to-use suspension.

For laboratory use only.

Applications

- Stop solution for TMB substrate
- ELISA

Specifications

Appearance	Clear liquid solution
Storage	4°C to 8°C
Shelf life	2 years

Directions for use

The product is a ready-to-use suspension. Use lab coat, gloves and goggles and handle with care since it contains acid.

Shipping and storage

The product is shipped at room temperature and is stable for 2 years when stored at 4-8°C. Store in original container, needs to be kept dark.

Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 and ISO 13485:2003 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



Ordering information

Article no.	Product name	Pack size
0F6EJ1 €€€ €	Stop solution for TMB	200 ml