





Manufactured By: HYPHEN BioMed

MICROPARTICLES

Novel assay for quantitatively measuring microparticles procoagulant activity in human citrated plasma. Microparticles are released from some blood cells, mainly platelet or endothelial cells, but also from monocytes or lymphocytes, upon cell activation. They frequently expose phosphatidyl serine, which focuses initiation of blood coagulation pathways. Microparticles are the markers arising as a consequence of the disease, but can also be the cause of worsening of the disease through their procoagulant activity.



Form AH59
5-2007



7768 Service Center Drive • West Chester OH 45069

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com