







# ANIARA

Manufactured By: U-CyTech Biosciences



## CYTOKINE ELISPOT & ELISA TECHNOLOGY



**ANIARA**

7768 Service Center Drive • West Chester OH 45069

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: [info@aniara.com](mailto:info@aniara.com)

[www.aniara.com](http://www.aniara.com)

Form AU01  
4-2007

# Cytokine ELISPOT and ELISA assay range

## ELISPOT kit; black spots

2-plate and 5-plate format

PVDF- or polystyrene-bottomed plates

- Coating antibodies
- Biotinylated detector antibodies
- $\phi$ -labeled anti-biotin antibodies
- Substrate
- Blocking stock solution
- Dilution buffer
- Tween-20
- With or without ELISPOT plates
- Detailed manual

## ELISPOT kit; red spots

2-plate, 5-plate and 20-plate format

PVDF-bottomed plates

- Coating antibodies
- Biotinylated detector antibodies
- Streptavidin-HRP conjugate
- Substrate
- Blocking stock solution
- Dilution buffer
- Tween-20
- With or without ELISPOT plates
- Detailed manual

## ELISA kit

5-plate format

- Coating antibodies
- Biotinylated detector antibodies
- Cytokine standards
- Streptavidin-HRP conjugate
- Substrate tablets + buffer capsules
- BSA stock solution
- Cytokine stabilization buffer
- Tween-20
- ELISA plates
- Detailed manual

Cytokine	ELISPOT range Species				ELISA range Species			
	Human	Monkey	Mouse	Rat	Human	Monkey	Mouse	Rat
GM-CSF	√	√			√	√		
Granzyme B	√	√						
IFN- $\gamma$	√	√	√	√	√	√	√	√
IL-1 $\beta$	√	√			√	√		
IL-2	√	√			√	√		
IL-4	√	√	√	√	√	√	√	√
IL-5	√	√	√		√	√		
IL-6	√	√			√	√		
IL-8	√	√						
IL-10	√	√	√		√	√		
IL-12p40+p70	√	√			√	√		
IL-12p70	√				√			
IL-13	√	√			√	√		
TNF- $\alpha$	√	√	√	√	√	√	√	√

## Also available

Matched antibody pairs for cytokine  
ELISPOT and ELISA analyses  
10- and 20-plate formats

## Guaranteed stability

We warrant that our antibody pairs for cytokine ELISPOT and ELISA analyses will remain stable for minimal one year from the date of purchase if stored properly. We will replace or refund the purchase price of any product returned to us that we find to be out of compliance with our specifications.

