

5D Protein Stabilizer

*For improved protein activity, shelf life,
stability, cell attachment and more.*



5-Diagnostics

Research and Diagnostic Solutions. Delivered.

ANIARA
DIAGNOSTICA

5-D | HYPHEN BIOMED | MEDICAGO | XENOMETRIX

Improve your protein & cell culture protocols with a BSA Alternative: 5-D Protein Stabilizer

5-D Protein Stabilizer, a porcine collagen peptide fraction, the alternative for BSA to improve your protein & cell culture protocols.

- Stabilizes and Increases the shelf life of proteins and enzymes
- Maintains biological activity of proteins during lyophilisation, freeze-thaw cycles, urea, or ethanol treatment
- Prevents protein adsorption on glass surfaces
- Stabilizes proteins during heat treatment
- BSA alternative with minimal lot to lot variation

Article No.	Product name	Quantity
A5D-82411A	5D Protein Stabilizer	100 mL
A5D-82411B	5D Protein Stabilizer	1000 mL

For laboratory use only



Local contact:

ANIARA
DIAGNOSTICA
5-D | HYPHEN BIOMED | MEDICAGO | XENOMETRIX

ANIARA Diagnostica
7768 Service Center Dr.
West Chester Township
OH, 45069, United States

Tel: +1 513-770-1991
info@aniara.com
www.aniara.com

5-Diagnosics Ltd.

Heuberg 7 | CH 4051 Basel | Switzerland | Tel: +41 61 551 551 4 | info@5diagnostics.com | www.5-diagnostics.com

***For Research Use Only. Not For Use In Diagnostic Procedures**