


Read carefully the instructions in the system user manual and on the labeling and instructions for use of the reagent. Manufactured by HYPHEN BioMed a Sysmex Group Company. Test details, information or availability varies according to country.
©2022 HYPHEN BioMed. All rights reserved. All trademarks are the property of HYPHEN BioMed unless otherwise specified. Sysmex is a trademark of Sysmex Corp. All other trademarks and brands are the property of their respective owners.
Local distributors' list on www.hyphen-biomed.com too.



ANIARA
DIAGNOSTICA

S-D | HYPHEN BIOMED | MEDICAGO | XENOMETRIX

7768 Service Center Dr.,
West Chester Township,
OH, 45069, United
States

Technical support:
info@aniara.com
www.aniara.com

H-002-00 - Photo Credit: Shutterstock



Direct Oral Anticoagulants (DOACs)

A COMPLETE LINE FOR ANTICOAGULANT DETECTION

ANIARA
DIAGNOSTICA

S-D | HYPHEN BIOMED | MEDICAGO | XENOMETRIX

Aniara Diagnostica is a distributor for Hyphen BioMed

Developed and optimized for your research laboratory applications

- Multiplatform
 - method application available.
- Time saving
 - High practicability and full automation,
 - Liquid presentation for the routine anti-FXa assay.
- High quality and reliable standards
 - Calibrators traceability to WHO International Standards, when available, or to Pharmacopoeias' reference preparations,
 - No interference of released heparin neutralizing proteins in the kinetics BIOPHEN™ Heparin LRT kit, thanks to the use of Dextran Sulfate,
 - Long shelf life of reagents.

Liquid, one-stage anti-FXa assay for heparins and FXa-DOACs

- Liquid, ready to use.
- Drug-dedicated plasma calibrators and controls.

REAGENT				CALIBRATORS		CONTROLS	
Product Name	RS	Ref.	Presentation	Product Name	Ref.	Product Name	Ref.
• BIOPHEN™ Heparin LRT Ready to use liquid reagents	RUO	A221011-RUO	R1: 4 x 7.5 mL	• BIOPHEN™ UFH Calibrator • BIOPHEN™ Heparin Calibrator • BIOPHEN™ Orgaran® Calibrator • BIOPHEN™ Arixtra® Calibrator • BIOPHEN™ Rivaroxaban Calibrator • BIOPHEN™ Apixaban Calibrator • BIOPHEN™ Edoxaban Calibrator	A222301	• BIOPHEN™ UFH Control • BIOPHEN™ LMWH Control	A223101
			R2: 4 x 7.5 mL		A2223901		A224101
					A222001		A223901
					A222201		A223701
					A222501-RUO		A224301
					A222701-RUO		A224401
		A221013-RUO	R1: 3 x 3 mL	• BIOPHEN™ Orgaran® Control • BIOPHEN™ Arixtra® Control Plasma • BIOPHEN™ Rivaroxaban Control • BIOPHEN™ Apixaban Control • BIOPHEN™ Edoxaban Control	A2226001-RUO	• BIOPHEN™ Orgaran® Control • BIOPHEN™ Arixtra® Control Plasma • BIOPHEN™ Rivaroxaban Control • BIOPHEN™ Apixaban Control • BIOPHEN™ Edoxaban Control	A223801
			R2: 3 x 3 mL		A226101-RUO		A224201
					A226201-RUO		A223501-RUO
					A226401-RUO		A224001-RUO
					A226501-RUO		A225101-RUO
					A225501-RUO		A224501-RUO
		A221015-RUO	R1: 4 x 5 mL				A225201-RUO
			R2: 4 x 5 mL				A225301-RUO
							A225401-RUO
							A225501-RUO

HIT/T Detection with ZYMUTEST HIA

- Dynamic immunoassay mimicking the “in vivo pathogenic mechanism”.
- ELISA method for the semi-quantitative detection of Heparin-dependent antibodies.
- Successfully reported in many publications for identification of anti-PF4 antibodies present in VITT.

Product Name	RS	Ref.	Presentation
• ZYMUTEST HIA IgG	CE-IVD 510(k)	ARK040A	96 tests

Product Name	RS	Ref.	Presentation
• ZYMUTEST HIA MonoStrip IgG	CE-IVD 510(k)	ARK041A	32 tests

Alternative DiXal dedicated assay for DOACs, insensitive to heparins

REAGENT				CALIBRATORS		CONTROLS	
Product Name	RS	Ref.	Presentation	Product Name	Ref.	Product Name	Ref.
• BIOPHEN™ DiXal	RUO	A221030-RUO	R1: 3 x 2.5 mL	• BIOPHEN™ Rivaroxaban Calibrator • BIOPHEN™ Apixaban Calibrator • BIOPHEN™ Edoxaban Calibrator	A222701-RUO A226001-RUO A226101-RUO A226201-RUO A226401-RUO A226501-RUO	• BIOPHEN™ Rivaroxaban Control • BIOPHEN™ Apixaban Control • BIOPHEN™ Edoxaban Control	A224501-RUO A225101-RUO A225201-RUO A225301-RUO A225401-RUO A225501-RUO
			R2: 3 x 2.5 mL				
			R3: 4 x 20 mL				

High Sensitive Anti-IIa Assays for Direct Thrombin Inhibitors

- Bivalirudin.
- Dabigatran.
- Argatroban.
- Drug-dedicated plasma calibrators and controls, traceable to International Standards / reference material.

REAGENT				CALIBRATORS		CONTROLS	
Product Name	RS	Ref.	Presentation	Product Name	Ref.	Product Name	Ref.
• BIOPHEN™ DTI	RUO	A220202-RUO	R1: 2 x 2.5 mL	• BIOPHEN™ Dabigatran Calibrator • BIOPHEN™ Bivalirudin Calibrator	A222901-RUO A222801-RUO A226701-RUO	• BIOPHEN™ Dabigatran Control • BIOPHEN™ Bivalirudin Control	A225001-RUO A224701-RUO A225701-RUO
			R2: 2 x 2.5 mL				
			R3: 2 x 25 mL				

REAGENT				CALIBRATORS		CONTROLS	
Product Name	RS	Ref.	Presentation	Product Name	Ref.	Product Name	Ref.
• HEMOCLOT™ Thrombin Inhibitors	RUO	ACK002K-RUO	R1: 3 x 1 mL	• BIOPHEN™ Dabigatran Calibrator • Argatroban Plasma Calibrator • BIOPHEN™ Bivalirudin Calibrator	A222901-RUO A222801-RUO ASC030K-RUO A226701-RUO	• BIOPHEN™ Dabigatran Control • Argatroban Control Plasma • BIOPHEN™ Bivalirudin Control	A225001-RUO A224701-RUO ASC035K-RUO A225701-RUO
			R2: 3 x 1 mL				
		ACK002L-RUO	R1: 3 x 2.5 mL				
			R2: 3 x 2.5 mL				



Auxiliary Tris-NaCl-EDTA buffer at pH 7.85 (ref. AAR032A-RUO-MIN, AAR032K-RUO) to allow measuring anti-FXa DOACs with BIOPHEN™ Heparin LRT kits (ref. A221011-RUO/A221013-RUO/A221015-RUO), without interference of heparin.