

— PLASMIN / PLASMINOGEN CHROMOGENIC SUBSTRATE · BIOPHEN CS · HYPHEN BIOMED

Measure plasmin / plasminogen — *read fibrinolysis directly.*

A highly purified synthetic peptide substrate for **plasmin / plasminogen**. On cleavage it releases **para-nitroaniline (pNA)**, read at 405 nm in direct proportion to enzyme activity. Built on the established **Chromogenix S-2403** chemistry and backed by HYPHEN BioMed's own head-to-head equivalence studies — supplied lyophilized, single vial to bulk. Part of the **BIOPHEN™ CS series**, distributed by Aniara in North America and the Nordic & Baltic states.

CHROMOGENIC SUBSTRATE · HYPHEN BIOMED

Plasmin Substrate

pNA at 405 nm · S-2403 equivalent

- Releases **para-nitroaniline (pNA)** on cleavage by plasmin / plasminogen — absorbance at 405 nm is proportional to activity.
- Peptide sequence **H-D-Val-Leu-Lys-pNA** — the S-2403 substrate, no protocol redesign.
- **BIOPHEN CS-41(03) (A229041)** — single-vial (1x25 mg) to bulk (12x25 mg).
- **Documented equivalence** to the legacy Chromogenix substrate — head-to-head data on file.

S-2403 equivalent · RUO

25 mg

DOCUMENTED EQUIVALENCE

HYPHEN BioMed has run **head-to-head characterization studies** of the CS series against the matching Chromogenix substrates — not just an equivalence claim.

Identical peptide chemistry and **pNA-at-405 nm** readout — drop into existing protocols with no redesign.

A reliable second source for labs and manufacturers built around the legacy Chromogenix S-series.

— ONE SUBSTRATE, THREE BENCHES

Academic research

Fibrinolysis research — plasmin and plasminogen (SK-activated) kinetics.

Quality control

Plasminogen-assay and plasmin-reagent QC.

Pharma research

Thrombolytic and plasminogen-activator potency; antifibrinolytic screening.

Equivalence data, on request

HYPHEN BioMed head-to-head comparison studies versus the matching Chromogenix substrate, plus peptide sequence, MW, K_m , k_{cat} and molarity per vial, are available for assay development and validation. The user establishes and validates assay conditions and protease specificity.

Regulatory & availability: Plasmin Substrate (BIOPHEN™ CS series) is a synthetic, highly purified **research-grade reagent** — for research and pharmaceutical quality-control applications, not for use in diagnostic procedures. The thrombin (anti-IIa) and Factor Xa (anti-Xa) substrates support heparin-potency QC under the European & US Pharmacopoeia. Manufactured by HYPHEN BioMed and distributed by Aniara Diagnostica in **North America and the Nordic & Baltic states**. “Chromogenix” and S-series designations are referenced for cross-identification only and remain the marks of their respective owner. Single-vial and bulk formats; pricing on request. © 2026 Aniara Diagnostica, Inc.