


## Typical data Rat TNF- $\alpha$ ELISA kit

Catalogue no.: CT075A (5-plate format)

### Validated for:

Rat



Example of a representative Rat TNF- $\alpha$  standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.  
Sensitivity: 2 pg/ml.  
Specificity: natural and recombinant rat TNF- $\alpha$ .  
Calibration: recombinant rat TNF- $\alpha$ .  
Type of sample: serum, plasma and culture supernatant.

### References list:

- Bakker, J.M., Kavelaars, A., Kamphuis, P.J., Cobelens, P.M., van Vugt, H.H., van Bel, F. and Heijnen, C.J. 2000. *J. Immunol.* 165: 5932-5937
- Dieleman, J.M., de Lange, F., Houston, R.J.F., Biessels, G-J., Bär, P.R., Mackensen, B.G., Grocott, H.P. and Kalkman, C.J. 2006. *Life Sci.* 79: 551-558
- Haerter, K., Vroon, A., Kavelaars, A., Heijnen, C.J., Limmroth, V., Espinosa, E., Schedlowski, M. and Elsenbruch, S. 2004. *J. Neuroimmunol.* 146: 126-132
- Nair, V., Singh, S. and Gupta, Y.K. 2013. *J. Ayurveda Intergr. Med.* 4: 13-18

More references can be found in our References Database on our website: [www.ucytech.com/references](http://www.ucytech.com/references).

PROUDLY DISTRIBUTED BY   
ANIARA DIAGNOSTICA LLC  
+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com