


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Monkey IFN- γ ELISPOT kit

Catalogue no.: CT121-PR2 (2-plate format)
CT121-PR5 (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (1×10^4 cells/well)

Example of IFN- γ specific spots produced by rhesus macaque PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

Zhao, L., Liu, B., Ren, J., Feng, J., Pang, Z., Gao, J., Zhang, H., Tan, W., Tian, H. and Ruan L. 2011. Virol. J. 15: 302

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com