

U-CyTech BV  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Monkey GM-CSF ELISPOT kit


Catalogue no.: CT124-PR2 (2-plate format)  
CT124-PR5 (5-plate format)

### Validated for:

Rhesus macaque, Cynomolgus monkey, Mandrill



Stimulus: none ( $2 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $2.5 \times 10^4$  cells/well)

Example of GM-CSF specific spots produced by cynomolgus monkey PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com