


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Monkey IL-13 ELISPOT kit


Catalogue no.: CT132-T2 (2-plate format)
CT132-T5 (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Japanese macaque, Crested black macaque, Lion-tailed macaque, Baboon, Mandrill, African green monkey, Black mangabey, Hanuman langur, Marmoset



Stimulus: none (5×10^4 cells/well)



Stimulus: concanavalin A (5×10^4 cells/well)

Example of IL-13 specific spots produced by rhesus macaque PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

Yoshino, N., Lü, F.X-S., Fujihashi, K., Hagiwara, Y., Kataoka, K., Lu, D., Hirst, L., Honda, M., van Ginkel, F.W., Takeda, Y., Miller, C.J., Kiyono, H. and McGhee, J.R. 2004. *J. Immunol.* 173: 6850-6857

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com