

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Monkey GM-CSF ELISA kit

Catalogue no.: CT140A (5-plate format)

Validated for:

Rhesus macaque and Cynomolgus monkey



Example of a representative Monkey GM-CSF standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.

Sensitivity: 2 pg/ml.

Calibration: recombinant and natural rhesus macaque GM-CSF.

Type of sample: serum, plasma and culture supernatant.

References list:

Montiel N.A., Todd, P.A., Yee, J.A. and Lerche N.W. 2012. Comp. Med. 62: 61-68

PROUDLY DISTRIBUTED BY 
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com