


Typical data Monkey IL-4 ELISA kit

Catalogue no.: CT143A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Barbary macaque, Baboon



Example of a representative Monkey IL-4 standard curve using a 4-parameter logistic fit.

Recommended assay range: 0.5-32 pg/ml.

Sensitivity: 0.5 pg/ml.

Calibration: natural rhesus macaque IL-4.

Type of sample: serum, plasma and culture supernatant.

No interference by the soluble IL-4 receptor at a 1000-fold molar excess.

References list:

- Fallon, P.G., Gibbons, J., Vervenne, R.A., Richardson, E.J., Fulford, A.J.C., Kiarie, S., Sturrock, R.F., Coulson, P.S., Deelder, A.M., Langermans, J.A.M., Thomas, A.W. and Dunne, D.W. 2003. *J. Infect. Dis.* 187: 939-945
- Hartmann, G., Marschner, A., Viveros, P.R., Stahl-Hennig, C., Eisenblätter, M., Suh, Y-S., Endres, S., Tenner-Racz, K., Überla, K., Racz, P., Steinman, R.M. and Ignatius, R. 2005. *Vaccine* 23: 3310-3317
- Heydenburg Fuller, D., Rajakumar, P.A., Wilson, L.A., Trichel, A.M., Fuller, J.T., Shipley, T., Wu, M.S., Weis, K., Rinaldo, C.R., Haynes, J.R. and Murphey-Corb, M. 2002. *J. Virol.* 76: 3309-3317
- Van Rompay, K.K.A., Marthas, M.L. and Bischofberger, N. 2004. *Antiviral Res.* 63: 133-138
- De Swart, R.L., Kuiken, T., Timmermans, H.H., van Amerongen, G., van den Hoogen, B.G., Vos H.W., Neijens, H.J., Andeweg, A.C. and Osterhaus, A.D.M.E. 2002. *J. Virol.* 76: 11561-11569
- Tenbusch, M., Ignatius, R., Nchinda, G., Trumfheller, C., Salazar, A.M., Töpfer, K., Saueremann, U., Wagner, R., Hannaman, D., Tenner-Racz, K., Racz, P., Stahl-Hennig, C. and Überla, K. 2012. *PLoS ONE* 7: e39038
- Zhang W., Ahmad G., Torben W., Noor Z., Le L., Damian R.T., Wolf R.F., White G.L., Chavez-Suarez M., Podesta R.B., Kennedy R.C., and Siddiqui A.A. 2010. *J. Infect. Dis.* 201: 1105-1112

More references can be found in our References Database on our website: www.ucytech.com/references.