


**U-CyTech BV**

Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

**Typical data Monkey IL-5 ELISA kit****Catalogue no.:** CT144A (5-plate format)**Validated for:**

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Barbary macaque, Baboon, African green monkey



Example of a representative Monkey IL-5 standard curve using a 4-parameter logistic fit.

Recommended assay range: 0.5-32 pg/ml.  
Sensitivity: 0.5 pg/ml.  
Calibration: natural rhesus macaque IL-5.  
Type of sample: serum, plasma and culture supernatant.

**References list:**

Fallon, P.G., Gibbons, J., Vervenne, R.A., Richardson, E.J., Fulford, A.J.C., Kiarie, S., Sturrock, R.F., Coulson, P.S., Deelder, A.M., Langermans, J.A.M., Thomas, A.W. and Dunne, D.W. 2003. J. Infect. Dis. 187: 939-945

Mascarell, L., Bauche, C., Fayolle, C., Diop, O.M., Dupuy, M., Nougarede, N., Perraut, R., Ladant, D. and Leclerc, C. 2006. Vaccine 24: 3490-3499

Polakos, N.K., Drane, D., Cox, J., Ng, P., Selby, M.J., Chien, D., O'Hagan, D.T., Houghton, M. and Paliard, X. 2001. J. Immunol. 166: 3589-3598

De Swart, R.L., Kuiken, T., Timmermans, H.H., van Amerongen, G., van den Hoogen, B.G., Vos H.W., Neijens, H.J., Andeweg, A.C. and Osterhaus, A.D.M.E. 2002. J.Virol. 76: 11561-11569

De Waal, L., Power, U.F., Yüksel, S., van Amerongen, G., Nguyen, T.N., Niesters, H.G.M., de Swart, R.L. and Osterhaus, A.D.M.E. 2004. Vaccine 22: 915-922

PROUDLY DISTRIBUTED BY   
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com