


## Typical data Monkey IL-6 ELISA kit

Catalogue no.: CT145A (5-plate format)

### Validated for:

Rhesus macaque, Cynomolgus monkey, Barbary macaque, Japanese macaque, Baboon, African green monkey



Example of a representative Monkey IL-6 standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.

Sensitivity: 2 pg/ml.

Calibration: natural rhesus macaque IL-6.

Type of sample: serum, plasma and culture supernatant.

Minor interference (<10%) by the sIL-6 receptor at a 1000-fold molar excess.

### References list:

Beignon, A.-S., Mollier, K., Liard, C., Coutant, F., Munier, S., Rivière, J., Souque, P. and Charneau, P. 2009. *J. Virol.* 83: 10963-10974

Coe, C.L., Lubach, G.R. and Kinnard, J. 2012. *Age (Dordr)* 34: 1169-1177

Jacquelin, B., Petitjean, G., Kunkel, D., Liovat, A.-S., Jochems, S.P., Rogers, K.A., Ploquin, M.J., Madec, Y., Barré-Sinoussi, F., Dereuddre-Bosquet, N., Lebon, P., Le Grand, R., Villinger, F. and Müller-Trutwin, M. 2014. *PLoS Pathog.* 10: e1004241

McFarlane, D., Wolf, R.F., McDaniel, K.A. and White, G.L. 2012. *J. Med. Primatol.* 41: 341-348

Rout, N., Greene, J., Yue, S., O'Connor, D., Johnson, P.R., Else, J.G., Exley, M.A. and Kaur, A. 2012. *PLoS Pathog.* 8: e1002928

More references can be found in our References Database on our website: [www.ucytech.com/references](http://www.ucytech.com/references).

PROUDLY DISTRIBUTED BY



+1 (513) 770-1991 | +1 (866) 783-3797  
 7768 Service Center Drive  
 West Chester, OH 45069  
[info@aniara.com](mailto:info@aniara.com)