


Typical data Monkey IL-12/23p40 ELISA kit

Catalogue no.: CT149A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey



Example of a representative Monkey IL-12/23p40 standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.

Sensitivity: 4 pg/ml.

Calibration: recombinant rhesus macaque IL-12p70.

Type of sample: serum, plasma and culture supernatant.

The assay detects free p40 monomers, p40 dimers and the p40 subunits of heterodimeric IL-12 (p35/p40) and heterodimeric IL-23 (p19/p40).

References list:

Brown, J.N., Ortiz, G.M., Angel, T.E., Jacobs, J.M., Gritsenko, M., Chan, E.Y., Purdy, D.E., Murnane, R.D., Larsen, K., Palermo, R.E., Shukia, A.K., Clauss, T.R., Katze, M.G., McCune, J.M. and Smith, R.D. 2012. Mol. Cell Proteomics 11: 605-618

Fallon, P.G., Gibbons, J., Vervenne, R.A., Richardson, E.J., Fulford, A.J.C., Kiarie, S., Sturrock, R.F., Coulson, P.S., Deelder, A.M., Langermans, J.A.M., Thomas, A.W. and Dunne, D.W. 2003. J. Infect. Dis. 187: 939-945

Jacquelin, B., Petitjean, G., Kunkel, D., Liovat, A.-S., Jochems, S.P., Rogers, K.A., Ploquin, M.J., Madec, Y., Barré-Sinoussi, F., Dereuddre-Bosquet, N., Lebon, P., Le Grand, R., Villinger, F. and Müller-Trutwin, M. 2014. PLoS Pathog. 10: e1004241

Mehlhop, E., Villamide, L.A., Frank, I., Gettie A., Santisteban, C., Messmer, D., Ignatius, R., Lifson, J.D. and Pope M. 2002. J. Immunol. Meth. 260: 219-234

Mortara, L., Ploquin, M. J.-Y., Faye, A., Scott-Algara, D., Vaslin, B., Butor, C., Hosmalin, A., Barré-Sinoussi, F., Diop, O. and Müller-Trutwin, M.C. 2006. J. Immunol. Meth. 308: 138-155

Van Rompay, K.K.A., Marthas, M.L. and Bischofberger, N. 2004. Antiviral Res. 63: 133-138

More references can be found in our References Database on our website: www.ucytech.com/references.