


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Monkey G-CSF ELISA kit

Catalogue no.: CT155A (5-plate format)

### Validated for:

Rhesus macaque



Example of a representative Monkey G-CSF standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.

Sensitivity: 2 pg/ml.

Calibration: recombinant and natural rhesus macaque G-CSF.

Type of sample: serum, plasma and culture supernatant.

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com