


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Monkey MCP-1 ELISA kit

Catalogue no.: CT156A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey



Example of a representative Monkey MCP-1 standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.

Sensitivity: 8 pg/ml.

Calibration: natural rhesus macaque MCP-1.

Type of sample: serum, plasma and culture supernatant.

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com