


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human GM-CSF ELISA kit

Catalogue no.: CT200A (5-plate format)

Validated for:

Human



Example of a representative Human GM-CSF standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.
Sensitivity: 2 pg/ml.
Specificity: natural and recombinant human GM-CSF.
Calibration: recombinant human GM-CSF.
Type of sample: serum, plasma and culture supernatant.

References list:

Duan, R., Remeijer, L., van Dun, J.M., Osterhaus, A.D.M.E. and Verjans, G.M.G.M. 2007. Invest. Ophthalmol. Vis. Sci. 48: 277-284

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com