


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-5 ELISA kit

Catalogue no.: CT204A (5-plate format)

Validated for:

Human



Example of a representative Human IL-5 standard curve using a 4-parameter logistic fit.

Recommended assay range: 0.5-32 pg/ml.
Sensitivity: 0.5 pg/ml.
Specificity: natural and recombinant human IL-5.
Calibration: recombinant human IL-5.
Type of sample: serum, plasma and culture supernatant.

References list:

Chiriva-Interati, M., Yu, Y., Mirandola, L., D'Cunha, N., Hardwicke, F., Cannon, M.J., Cobos, E. and Kast M.W. 2012. Prostate 72: 12-23

Eberl, M., Langermans, J.A., Frost, P.A., Vervenne, R.A., van Dam, G.J., Deelder, A.M., Thomas, A.W., Coulson, P.S. and Wilson, R.A. 2001. Infect. Immun. 69: 5352-5362

More references can be found in our References Database on our website: www.ucytech.com/references.

PROUDLY DISTRIBUTED BY



ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com