


**U-CyTech BV**

Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

**Typical data Human TNF- $\alpha$  ELISA kit****Catalogue no.:** CT209A (5-plate format)**Validated for:**

Human



Example of a representative Human TNF- $\alpha$  standard curve using a 4-parameter logistic fit.

Recommended assay range: 1-64 pg/ml.  
Sensitivity: 1 pg/ml.  
Specificity: natural and recombinant human TNF- $\alpha$ .  
Calibration: recombinant human TNF- $\alpha$ .  
Type of sample: serum, plasma and culture supernatant.  
No interference by the sTNF receptors (RI/RII).

**References list:**

- Chang, W.L.W., Baumgarth, N., Yu, D. and Barry, P.A. 2004. J. Virol. 78: 8720-8731
- Chiriva-Interati, M., Mirandola, L., D'Cunha, N., Hardwicke, F., Cannon, M.J., Cobos, E. and Kast, M.W. 2012. Prostate 72: 12-23
- Jagessar S.A., Kap Y.S., Heijmans N., van Driel N., van Straalen L., Bajramovic J.J., Brok H.P., Blezer E.L., Bauer J., Laman J.D. and 't Hart B.A. 2010. J. Neuropathol. Exp. Neurol. 69: 372-385
- Skov, L., Allen, M.H., Bang, B., Francis, D., Barker, J.N.W.N. and Baadsgaard, O. 2003. Exp. Dermatol. 12: 772-776
- Tang, Y-M., Chen, M-H., Chen, G-H., Cai, C-J., He, X-S., Lu, M-G. and Bao, W-M. 2005. World J. Gastroenterol. 11: 4574-4578
- Weinhold, M., Eisenblätter, M., Jasny, E., Fehlings, M., Finke, A., Gayum, H., Rüschen-dorf, U., Renner Viveros, P., Moos, V., Allers, K., Schneider, T., Schaible, U.E., Schumann, R.R., Mielke, M.E. and Ignatius, R. 2013. PLoS ONE 8: e65934

More references can be found in our References Database on our website: [www.ucytech.com/references](http://www.ucytech.com/references).

PROUDLY DISTRIBUTED BY   
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
[info@aniara.com](mailto:info@aniara.com)