

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Human IL-12p70 ELISA kit

Catalogue no.: CT210A (5-plate format)

Validated for:

Human



Example of a representative Human IL-12p70 standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.

Sensitivity: 4 pg/ml.

Specificity: natural and recombinant human IL-12.

Calibration: recombinant human IL-12p70.

Type of sample: serum, plasma and culture supernatant.

The assay detects human IL-12p70 and not free p40 or p35 subunits.

References list:

Hoogendoorn, M., Olde Wolbers, J., Smit, W.M., Schaafsma, M.R., Jedema, I., Barge, R.M.Y., Willemze, R. and Falkenburg, J.H.F. 2005. Clin. Cancer Res. 11: 5310-5318

Nauta, A.J., Krusselbrink, A.B., Lurvink, E., Willemze, R. and Fibbe, W.E. 2006. J. Immunol. 177: 2080-2087

Verkade, M.A., van Druningen, C.J., Vaessen, L.M.B., Hesselink, D.A., Weimar, W. and Betjes, M.G.H. 2007. Nephrol. Dial. Transplant. 22: 128-138

Wang, H., Su, X., Zhang, P., Liang, J., Wei, H., Wan, M., Wu, X., Yu, Y. and Wang, L. 2011. Vaccine 29: 2328-2335

Weinhold, M., Eisenblätter, M., Jasny, E., Fehlings, M., Finke, A., Gayum, H., Rüschemdorf, U., Renner Viveros, P., Moos, V., Allers, K., Schneider, T., Schaible, U.E., Schumann, R.R., Mielke, M.E. and Ignatius, R. 2013. Plos ONE 8: e65934

More references can be found in our References Database on our website: www.ucytech.com/references.

PROUDLY DISTRIBUTED BY 
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com