


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Human IFN- $\gamma$ ELISPOT kit


Catalogue no.: CT230-PB2 (2-plate format)  
CT230-PB5 (5-plate format)

### Validated for:

Human



Stimulus: none ( $4 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $1 \times 10^4$  cells/well)

Example of IFN- $\gamma$  specific spots produced by human PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.

### References list:

Chiriva-Interati, M., Mirandola, L., Yu, Y., Jenkins, M.R., Gornati, R., Bernardini, G., Gioia, M, Chiamonte, R., Cannon, M.J., Kast, W.M. and Cobos, E. 2011. *J. Immunother.* 34: 490-499

Correale, P., Botta, C., Cusi, M., Del Vecchio, M., De Santi, M., Gori Savellini, G., Bestoso, E., Apollinari, S., Mannucci, S., Marra, M., Abruzzese, A., Aquino, A., Turriziani, M., Bonmassar, L., Caraglia, M. and Tagliaferri, P. 2011. *Int. J. Cancer Epub*

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com