


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-4 ELISPOT kit


Catalogue no.: CT232-PR2 (2-plate format)
CT232-PR5 (5-plate format)

Validated for:

Human



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (2.5×10^4 cells/well)

Example of IL-4 specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

Xue, X-K., Sun, W-J., Hu, J., Yao, H-P. and Cai, Z. 2006. Chin. Med. J. 119: 1038-1041