


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Human IL-5 ELISPOT kit


Catalogue no.: CT233-T2 (2-plate format)  
CT233-T5 (5-plate format)

### Validated for:

Human



Stimulus: none ( $2 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $2 \times 10^5$  cells/well)

Example of IL-5 specific spots produced by human PBMC.  
Preincubation: 42h, incubation, ELISPOT: 5h.

### References list:

Martin, S., Wolf-Eichbaum, D., Duinkerken, G., Scherbaum, W., Kolb, H., Noordzij J.G. and Roep, B.O. 2001. N. Engl. J. Med. 345: 1036-1040

Pfleger, C., Meierhoff, G., Kolb, H. and Schloot, N.C. 2010. J. Autoimmun. 34: 127-135

Schloot, N.C., Meierhoff, G., Faresjö, M.K., Ott, P., Putnam, A., Lehmann, P., Gottlieb, P., Roep, B.O., Peakman, M., Tree, T. 2003. J. Autoimmun. 21: 365-376

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com