


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-13 ELISPOT kit


Catalogue no.: CT236-PR2 (2-plate format)
CT236-PR5 (5-plate format)

Validated for:

Human



Stimulus: none (4×10^4 cells/well)



Stimulus: PMA/ionomycin (2.5×10^4 cells/well)

Example of IL-13 specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 140220