


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human TNF- α ELISPOT kit


Catalogue no.: CT237-PR2 (2-plate format)
CT237-PR5 (5-plate format)

Validated for:

Human



Stimulus: none (5×10^3 cells/well)



Stimulus: PMA/ionomycin (5×10^3 cells/well)

Example of TNF- α specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com