


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human GM-CSF ELISPOT kit


Catalogue no.: CT241-PB2 (2-plate format)
CT241-PB5 (5-plate format)

Validated for:

Human



Stimulus: none (2×10^5 cells/well)



Stimulus: PMA/ionomycin (1×10^4 cells/well)

Example of GM-CSF specific spots produced by human PBMC.
Incubation ELISPOT: 20h (no preincubation).



PROUDLY DISTRIBUTED BY **ANIARA**
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com