


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Human IL-1 $\beta$ ELISPOT kit


**Catalogue no.:** CT242-PR2 (2-plate format)  
CT242-PR5 (5-plate format)

### Validated for:

Human



Stimulus: none ( $5 \times 10^4$  cells/well)



Stimulus: PMA/ionomycin ( $5 \times 10^4$  cells/well)

Example of IL-1 $\beta$  specific spots produced by human PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com