

U-CyTech BV  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Data sheet    Monoclonal antibody to human interleukin 5 (IL-5)

Cat. No.:	CT255
Clone:	QS-5
Isotype:	Mouse IgG <sub>1</sub>
Production:	<i>In vitro</i> using serum free medium.
Purification:	Ion exchange chromatography.
Packaging:	Lyophilized and vacuum-packed.
Contents:	0.5 mg/vial
Buffer:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Specificity:	Binds with high efficiency to natural and recombinant human IL-5.
Sterility:	Membrane filtered (0.2 µm).
Reconstitution:	Dissolve the contents of the vial by injection of 0.5 ml sterile distilled water.
Stability:	Lyophilized product is stable for at least one year at 4°C (expiry date is indicated on the vial). After reconstitution, the contents can be safely stored at 4°C for one month or for one year at -20°C. Add 0.02% sodiumazide to prevent bacterial growth.
Quantitation:	Antibody concentration was determined by absorbance, taking A280=1.4 for a 1 mg/ml solution.
Application:	Coating antibody in ELISA Intracellular staining <i>In vitro</i> neutralization Western blot analysis



PROUDLY DISTRIBUTED BY



ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com