

U-CyTech BV


Yalelaan 48
 3584 CM Utrecht
 The Netherlands
 P +31.30.253 5960
 F +31.30.253 9344
 INFO@ucytech.com
 www.ucytech.com

Typical data Mouse IFN- γ ELISA kit

Catalogue no.: CT301A (5-plate format)

Validated for:

Mouse



Example of a representative Mouse IFN- γ standard curve using a 4-parameter logistic fit.

Recommended assay range: 1-64 pg/ml.
 Sensitivity: 1 pg/ml.
 Specificity: natural and recombinant mouse IFN- γ .
 Calibration: recombinant mouse IFN- γ .
 Type of sample: serum, plasma and culture supernatant.

References list:

Eijkelkamp, N., Cobelens, P.M., Sanders, V.M., Heijnen, C.J. and Kavelaars, A. 2004. J. Neuroimmun. 150: 3-9

Niemi, M-B., Pacheco-López, G., Kou, W., Härtling, M., del Rey, A., Besedovsky, H.O. and Schedlowski, M. 2006. Brain Behav. Immun. 20: 527-531

Pouw N., Treffers-Westerlaken E., Kraan J., Wittink F., Ten Hagen T., Verweij J. and Debets R. 2010. Cancer Immunol. Immunother. 59: 921-931

Vroon A., Kavelaars, A., Limmroth, V., Lombardi, M.S., Goebel, M.U., van Dam, A-M., Caron, M.G., Schedlowski, M. and Heijnen, C.J. 2005. J. Immunol. 174: 4400-4406

Zhao, Y., Liu, M.Y., Wang, X.L., Liu, X.L., Yang, Y., Zou, H.B., Sun, S.M., Yu, L., Rosenthal, B., Shi, H.N., Boireau, P. and Wu, X.P. 2013. Vet Parasitol. 194: 211-216

More references can be found in our References Database on our website: www.ucytech.com/references.

PROUDLY DISTRIBUTED BY 
 ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
 7768 Service Center Drive
 West Chester, OH 45069
 info@aniara.com