

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Mouse IL-10 ELISA kit

Catalogue no.: CT307A (5-plate format)

Validated for:

Mouse



Example of a representative Mouse IL-10 standard curve using a 4-parameter logistic fit.

Recommended assay range: 32-2048 pg/ml.

Sensitivity: 32 pg/ml.

Specificity: natural and recombinant mouse IL-10.

Calibration: recombinant mouse IL-10.

Type of sample: serum, plasma and culture supernatant.

References list:

Ghotloo, S., Haji Molla Hoseini, M., Alimohammadian, M.H., Khaze, V., Memarnejadian, A. and Rostami, A. 2015. Parasitol. Int. 64: 219-221

Hong, C., Zhang, T. and Gao, X-M. 2013. Mol. Immunol. 55: 237-246

Pouw, N., Treffers-Westertaken E., Kraan J., Wittink F., Ten Hagen T., Verweij J. and Debets R. 2010. Cancer Immunol. Immunother. 59: 921-931

Willemen, H.L.D.M., Eijkelkamp, N., Garza Carbajal, A., Wang, H., Mack, M., Zijlstra, J., Heijnen, C.J. and Kavelaars, A. 2014. J. Pain 15: 496-506

Zhao, Y., Liu, M.Y., Wang, X.L., Liu, X.L., Yang, Y., Zou, H.B., Sun, S.M., Yu, L., Rosenthal, B., Shi, H.N., Boireau, P. and Wu, X.P. 2013. Vet. Parasitol. 194: 211-216

More references can be found in our References Database on our website: www.ucytech.com/references.

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797

7768 Service Center Drive

West Chester, OH 45069

info@aniara.com

Version 160708