


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Mouse IL-2 ELISA kit

Catalogue no.: CT309A (5-plate format)

### Validated for:

Mouse



Example of a representative Mouse IL-2 standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.  
Sensitivity: 4 pg/ml.  
Specificity: natural and recombinant mouse IL-2.  
Calibration: recombinant mouse IL-2.  
Type of sample: serum, plasma and culture supernatant.

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com

Version 140822