


Typical data Mouse IFN- γ ELISPOT kit


Catalogue no.: CT317-PR2 (2-plate format)
CT317-PR5 (5-plate format)

Validated for:

Mouse



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (2.5×10^4 cells/well)

Example of IFN- γ specific spots produced by mouse (Balb/c) splenocytes.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

- Bai, B., Lu, X., Meng, J., Hu, Q., Mao, P., Lu, B., Chen, Z., Yuan, Z. and Wang, H. 2008. *Mol. Immunol.* 45: 868-875
- Hu, H., Lu, X., Tao, L., Bai, B., Zhang, Z., Chen, Y., Zheng, F., Chen, J., Chen, Z. and Wang, H. 2007. *Clin. Vaccine Immunol.* 14: 894-901
- Li, S., Peng, L., Zhao, W., Zhong, H., Zhang, F., Yan, Z. and Cao, H. 2011. *Vaccine* 29: 3695-3702
- Xu, Y., Wang, B., Chen, J., Wang, Q., Zhu, B., Shen, H., Qie, Y., Wang, J. and Wang, H. 2006. *Scand. J. Immunol.* 64: 476-481
- Zhang, Y., Zan, Y., Shan, M., Liu, C., Shi, M., Li, W., Zhang, Z., Liu, N., Wang, F., Zhong, W., Liao, F., Gao, G.F. and Tien, P. 2006. *Biochem. Biophys. Res. Commun.* 344: 581-587

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 110810