


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Mouse IL-4 ELISPOT kit


Catalogue no.: CT319-T2 (2-plate format)
CT319-T5 (5-plate format)

Validated for:

Mouse



Stimulus: none (2×10^4 cells/well)



Stimulus: Concanavalin A (2×10^4 cells/well)

Example of IL-4 specific spots produced by mouse (Balb/c) splenocytes.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

- Cai, Q., Peng, G., Bu, L., Lin, Y., Zhang, L., Lustigmen, S. and Wang 2007. *Vaccine* 25: 5155-5165
- Han, G., Li, Y., Wang, J., Wang, R., Chen, G., Song, L., Xu, R., Yu, M., Wu, Z., Qian, J. and Shen, B. 2005. *J. Immunol.* 174: 4516-4524
- Han, G., Wang, R., Chen, G., Wang, J., Xu, R., Feng, J., Yu, M., Wu, X., Qian, J., Shen, B. and Li, Y. 2008. *Immunol. Lett.* 115: 110-116
- Li, P., Wan, Q., Feng, Y., Liu, M., Wu, J., Chen, X. and Zhang, X-L. 2007. *Vaccine* 25: 1544-1551
- Wu, F., Yuan, X.Y., Huang, W.S. and Chen, Y.H. 2009. *Vaccine* 27: 6095-6101
- Yu, H., Jiang, L-F., Fang, D-Y., Yan, H-J., Zhou, J-J., Zhou, J-M., Liang, Y., Gao, Y., Zhao, W. and Long, B-G. 2007. *Virology* 359: 264-274

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 110530