


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Mouse IL-10 ELISPOT kit


Catalogue no.: CT320-PB2 (2-plate format)
CT320-PB5 (5-plate format)

Validated for:

Mouse



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (2×10^5 cells/well)

Example of IL-10 specific spots produced by mouse (CD1) splenocytes.
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 110810