


**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Mouse IL-5 ELISPOT kit


Catalogue no.: CT321-PB2 (2-plate format)  
CT321-PB5 (5-plate format)

### Validated for:

Mouse



Stimulus: none ( $4 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $2 \times 10^5$  cells/well)

Example of IL-5 specific spots produced by mouse (Balb/c) splenocytes.  
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com