


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Mouse TNF- α ELISPOT kit


Catalogue no.: CT322-T2 (2-plate format)
CT322-T5 (5-plate format)

Validated for:

Mouse



Stimulus: none (1×10^5 cells/well)



Stimulus: Concanavalin A (5×10^4 cells/well)

Example of TNF- α specific spots produced by mouse (CD-1) splenocytes.
Preincubation: 24h, incubation ELISPOT: 20h.

References list:

Chiriva-Internati, M., Yu, Y., Mirandola, L., Jenkins, M.R., Chapman, C., Cannon, M., Cobos, E. and Kast, W.M. 2010. PLoS One 5: e10471

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com