

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Marmoset TNF- α ELISA kit

Catalogue no.: CT342A (5-plate format)

Validated for:

Marmoset



Example of a representative Marmoset TNF- α standard curve using a 4-parameter logistic fit.

Recommended assay range: 1-64 pg/ml.

Sensitivity: 1 pg/ml.

Calibration: natural marmoset TNF- α .

Type of sample: serum, plasma and culture supernatant.

No interference by the sTNF receptors (RI/RII) at a 1000-fold molar excess.

References list:

Jagessar, A.S., Heijmans, N., Bauer, J., Blezer, E.L.A., Laman, J.D., Migone, T.-S., Devalaraja, M.N. and 't Hart, B.A. 2012. J. Neuroimmune Pharmacol. 7: 557-570

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 140820