

Data sheet 96-well Transparent polystyrene ELISA plates

Cat. No.:	CT361
Specificity:	These 96-well transparent polystyrene F-bottomed plates (Greiner Bio-One GmbH) have been selected for use in U-CyTech ELISA procedures.
Application:	The plates do not require a pre-treatment and the coating antibodies can directly be applied to the wells.
Contents:	6 plates with 10 adhesive cover slips.
Sterility:	-
Packaging:	The plates are packed per 2.
Storage:	Stable until the expiry date (indicated on the label) when stored at room temperature.
General:	Transparent colorless polystyrene plate.
Safety:	The product is not classified as hazardous according to Regulation (EC) no. 1272/2008 and to Directive 1999/45/EC and its amendments.

